


Application of Collision/Reaction Cell ICP-MS to the Analysis of Variable and Unknown Samples Without Requiring Matrix-Specific or Element-Specific Tuning and Cell Conditions

Introduction. ICP-MS is a multi-element analytical technique, which uses a high temperature plasma source to decompose, dissociate, atomise and ionise the elements of interest in a sample, providing high sensitivity, rapid analysis, low detection limits, simple spectra and routine operation. In simple solutions, few interferences obscure the measurement of the trace elements of interest; primarily the affected elements are those that are overlapped by polyatomic ions derived from the sample solvent and the plasma gas (O, H, Ar, C, N). In sample matrices of the types commonly analysed in routine laboratories, the situation with respect to interferences can be much more complicated, especially if, as is typically the case, the matrix is complex, variable and of unknown composition. Under these circumstances, it is essential that the ICP-MS is able to handle the variable sample matrices, with the minimum of information on the sample composition, in order to permit unknown samples to be analysed routinely.

Conventional ICP-MS

Numerous reports have shown that conventional ICP-MS (without collision/reaction cell) can be optimised to provide dramatic improvement in the capacity to measure "difficult" trace elements in the presence of high levels of a sample matrix, simply through the use of design and operating conditions which lead to a high and stable plasma temperature to ensure efficient sample decomposition and analyte ionisation. The calibration on the right shows As at single ng/mL levels in a matrix of 5.6% aqua regia, used for soil digestion.




Calibration for As (2 - 40ng/mL) in 5.6% aqua regia (4.2% HCl, 1.4% HNO₃) - Agilent 7500a (non-cell instrument)


Collision/Reaction Cell ICP-MS

The collision/reaction cell provides the analyst with an additional method for reducing or removing interfering species from the mass spectrum. In **reaction mode** (shown here), specific target interferences are removed from specific analyte masses, using well-defined reaction pathways. **Collision mode** (using pure He cell gas) is non-specific: it attenuates ALL polyatomic species using a single set of cell conditions.

In collision mode, the larger polyatomic (interferant) species undergo more collisions (and so lose more energy) than the smaller, monatomic (analyte) species. Interference removal occurs by energy discrimination and is therefore non-specific, but relies on a very small ion energy spread of all ions entering the cell



H₂ cell gas optimisation plot for ⁸⁰Se



7500cs Octopole Reaction Cell Schematic

Reported limitations of Cell-based ICP-MS Utilising Highly Reactive Gases Include:

- Loss of sensitivity (due to reaction) for some required analyte elements; may require several sets of gas conditions for full analysis, so increasing method development
- Inability of a single reaction gas to remove multiple interfering species on a single analyte (e.g. ⁸⁰Se in a complex matrix).
- Appearance of new interfering species (eg cluster ions); requires each matrix to be evaluated individually for new interfering species, before cell conditions for each analyte isotope can be determined – makes cell conditions matrix specific
- Multiple anomalous peaks due to traces of residual gas/matrix contamination in cell (providing reactive analyte loss and cluster ion formation, even in non-gas mode)
- Residual or new (cell-formed) interferences, leading to a requirement to utilise multiple, matrix-dependent interference correction equations
- Cell-formed interferences can interfere with internal standard - incorrect quantification

Due to the above limitations, the operating conditions required for accurate analysis of complex, high matrix samples are highly analyte and matrix dependent, effectively preventing the reliable analysis of unknown sample types in cell mode with a reactive gas.

Cell Operation with an Inert Cell Gas

Using an inert cell gas and efficient energy discrimination (ED) removes all of these reported problems, since the cell is operating as a non-reactive environment. The use of an inert gas allows ED to be used as the primary interference removal process, which eliminates the method development inherent in the use of reactive cell gases. This means that the analysis conditions are independent of the analyte being measured (and of the interferences being removed), multi-element analysis is possible and the operating conditions do not have to be optimised for each matrix type.

A further benefit of collision mode is that no interference correction equations are necessary. Use of ED depends on the use of a high purity inert cell gas, which is free from organic contaminants and water vapour, so plumbing of the gas supply requires care. A light collision gas (He) is preferred, since less analyte ion scattering occurs.

Ed McCurdy, Don Potter and Glenn Woods,




Agilent Technologies UK Ltd, Lakeside, Cheadle Royal Business Park, Stockport, Cheshire, SK8 3GR UK
Steven Wilbur, Agilent Technologies Inc, 3380 146th Place SE Suite 300, Bellevue, WA 98007, USA.

Mode of Operation

The operation of a cell in collision mode, using ED as the method of interference removal, requires strict control of ion energy entering the cell. This is achieved using a shielded torch configuration, which provides much smaller ion energy spread than end-grounded or centre-grounded unshielded torch configurations. The control of ion energy is essential for effective interference removal by ED, since it is the small difference in ion energy that results from the larger number of collisions that a polyatomic ion suffers, compared to a monatomic ion, which is used to distinguish between the analyte and interferant species at the cell exit. For this discrimination to be possible, the analyte and interferant ions must enter the cell with essentially the same energy.

Multi-element Analysis in Single Component Matrix

The spectra below show the ability of collision mode to remove multiple polyatomic ion species in a matrix of 2% H₂SO₄. In addition to the normal ArO and Ar₂ background species, SO, S₂, SO₂, ArS, S₂O and SO₃ interferences interfere with the 1st row transition elements, As and Se. Efficient removal of all S-based species is not possible in reaction mode. One reported approach for reaction mode proposed to measure analytes indirectly via the formation of cluster ions with NH₃ reaction gas. This approach is inherently unreliable in unknown complex matrices: He collision mode is with ED is the only reliable method for the measurement of Zn in a sulphate-containing matrix.




While 2% H₂SO₄ is a single-component matrix, this illustrates that He collision mode with ED removes multiple interferences using a single set of operating conditions, demonstrating its applicability to unknown, complex samples.

Capacity to Remove Multiple Interferences from Multiple Analytes

It is the ion filtering effect of the inert cell gas, combined with a small ion energy spread, which enables the removal of ANY polyatomic ion from the ion beam. This is because ALL polyatomic ions have substantially larger ionic cross-sections than a monatomic (analyte) ions.

Since the number of collisions that an ion suffers will be dependent on its ionic cross-section (assuming equal ion energy at the cell entrance), ALL polyatomic ions will collide more frequently than monatomic analyte ions and so will suffer greater collisional deceleration, allowing them to be rejected from the ion beam through the use of a positive voltage step (ED) at the cell exit. The identity of the polyatomic ion is irrelevant.



Analytes can have multiple interferences in complex sample matrices

Reliable Semi-quantitative Analysis is dependent on the maintenance of a predictable relationship between counts vs concentration for each analyte isotope. Using a simple, non-specific calibration based on the stored or measured mass/response curve for the instrument and using only constant corrections for isotopic abundance and degree of ionisation, semi-quantitative analysis can produce usable concentrations for many elements, based on calibrations that contain only a few reference elements

Element	NIST 1640 Traces in Water		SLRS-4 River Water		LGC 6010 Hard Drinking Water		LGC 6177 Landfill Leachate	
	Mass	SQ Conc.	Ref	Mass	SQ Conc.	Ref	Mass	SQ Conc.
Li 7	7	50.7	N/A	7	50.7	N/A	7	50.7
Be 9	9	35	34.94	9	0.007	<0.2200	9	<2.400
B 10	10	280	301.1	10	6	N/A	10	6.700
Na 23	23	30,000	29,350	23	2,400	20,000	23	1,500,000
Mg 24	24	5,700	5,819	24	1,600	3,700	24	62,000
Al 27	27	54	52	27	66	160	27	110
Si 28	28	4,800	4,730	28	1,900	4,300	28	22,000
P 31	31	35	N/A	31	41	670	31	11,600
S 34	34	110	N/A	34	110	470	34	1,600
K 39	39	1,000	994	39	680	5,100	39	810,000
Ca 44	44	9,300	7,045	44	5,900	6,200	44	77,000
Sc 45	45	<5.800E-3	N/A	45	<6.800E-3	0.0065	45	0.21
Ti 49	49	0.089	N/A	49	0.1	0.1	49	18
V 51	51	13	12.90	51	0.28	0.92	51	63
Cr 52	52	37	36.6	52	5	51	48	160
Mn 55	55	120	121.5	55	3.5	3.37	55	130
Fe 56	56	29	34.3	56	120	240	56	3,300
Co 59	59	19	20.28	59	0.045	0.033	59	20
Ni 60	60	26	27.4	60	1	0.67	60	170
Cu 63	63	87	86.2	63	2.2	1.81	63	41
Zn 66	66	86	83.2	66	1.3	540	66	260
Ge 72	72	32	N/A	72	2.8	N/A	72	130
As 75	75	24	26.67	75	0.1	0.68	75	55
Se 78	78	21	21.96	78	<1.400	N/A	78	13
Br 79	79	<0.1600	N/A	79	8.9	N/A	79	6,900
Rb 85	85	1.8	2	85	1.2	N/A	85	400
Y 89	89	120	124.2	89	26	26.3	89	980
Zr 90	90	0.091	N/A	90	0.047	N/A	90	49
Nb 93	93	0.0035	N/A	93	<1.200E-3	N/A	93	1.4
Mo 95	95	16	45.75	95	0.3	0.21	95	6.6
Ru 101	101	<3.400E-3	N/A	101	<3.500E-3	N/A	101	<0.04200
Rh 103	103	IS Element	N/A	103	IS Element	N/A	103	IS Element
Pd 106	106	<2.900E-3	N/A	106	<2.900E-3	N/A	106	<0.03500
Ag 107	107	7.7	7.62	107	0.011	N/A	107	4.3
Cd 114	114	22	22.79	114	0.011	0.012	114	0.94
In 115	115	0.0046	N/A	115	<1.500E-3	N/A	115	0.18
Sb 121	121	2.1	N/A	121	0.042	N/A	121	48
Te 125	125	15	13.79	125	0.32	0.23	125	11.9
I 127	127	<0.1200	N/A	127	<0.1200	N/A	127	3
Cs 133	133	0.078	N/A	133	0.014	N/A	133	1,200
Ba 137	137	140	140	137	12	110	137	3.5
La 139	139	0.42	N/A	139	0.35	N/A	139	0.24
Ce 140	140	0.52	N/A	140	0.48	N/A	140	0.7
Pr 141	141	0.076	N/A	141	0.062	N/A	141	<1.000E-2
Nd 146	146	0.35	N/A	146	0.15	N/A	146	0.2
Sm 147	147	0.072	N/A	147	0.027	N/A	147	<0.05000
Eu 153	153	0.0098	N/A	153	0.0065	N/A	153	0.036
Gd 157	157	<1.600E-3	N/A	157	0.032	N/A	157	0.036
Tb 159	159	0.0017	N/A	159	0.0017	N/A	159	0.0098
Dy 163	163	0.021	N/A	163	0.013	N/A	163	0.14
Ho 165	165	0.006	N/A	165	0.0024	N/A	165	0.061
Er 168	168	0.016	N/A	168	0.0084	N/A	168	0.052
Tm 169	169	<3.600E-4	N/A	169	0.001	N/A	169	0.011
Yb 172	172	0.0082	N/A	172	0.0054	N/A	172	0.044
Lu 175	175	6.50E-04	N/A	175	3.20E-04	N/A	175	0.01
Hf 178	178	<1.400E-3	N/A	178	0.0014	N/A	178	0.44
Ta 181	181	3.10E-04	N/A	181	3.10E-04	N/A	181	0.0066
W 182	182	0.017	N/A	182	0.0058	N/A	182	58
Re 185	185	0.0087	N/A	185	0.011	N/A	185	0.061
Os 190	190	<1.600E-3	N/A	190	0.0033	N/A	190	0.0034
Ir 193	193	IS Element	N/A	193	IS Element	N/A	193	IS Element
Pt 195	195	<1.300E-3	N/A	195	<1.300E-3	N/A	195	<0.01300
Au 197	197	0.0066	N/A	197	<1.000E-3	N/A	197	5.4
Hg 202	202	0.012	N/A	202	0.065	N/A	202	0.28
Tl 205	205	0.035	N/A	205	0.045	N/A	205	<7.600E-3
Pb 208	208	27	27.89	208	0.16	0.086	208	17
Bi 209	209	0.0015	N/A	209	0.001	N/A	209	2.4
Th 232	232	0.16	N/A	232	0.065	N/A	232	0.1
U 238	238	0.85	N/A	238	0.09	0.05	238	0.35

Analysis using Agilent 7500c in He Mode
All Elements under same analysis conditions
All Elements in ug/L

The table above shows the accurate determination of multiple elements in a variety of reference material matrices, using only a single set of He cell gas conditions for ALL elements in ALL samples. No interference correction equations were used.